OMRON

OMRON Automatic Blood Pressure Monitor HBP-9030



"Accurate Measurement Support" features to provide accurate measurements

HBP-9030 is equipped with "Accurate Measurement Support" features helping users easily take measurements by themselves.

"Elbow Detection Sensor" detects an elbow and indicates the correct arm position.

HBP-9030 adopts "Movable Cuff" which adjusts the best angle for all users and maintain the correct sitting posture naturally. These features allow more accurate measurements.





OMRON Automatic Blood Pressure Monitor HBP-9030

Easily replaceable Arm Unit



Replaceable Arm Unit is prepared as an optional accessary. Users can quickly replace Arm Unit in the installation place. HBP-9030 shows the timing to replace Arm Unit.



Arm Unit with a cuff cover



"Arm Unit icon" notifies the timing to replace.



"Accurate Measurement Support" guiding users to the correct sitting posture

Accurate Measurement Navigation

The Illustration and audio guide support accurate measurements.



© Elbow Detection Sensor

Guide user's arm to the correct position for measurement.





Movable cuff

The cuff provides angle control that supports taller people to have correct sitting and arm position comfortably.



USB Data Transfer

Measurement data can be directly transmitted via USB by simply connecting HBP-9030 to the connector of a terminal equipment such as a computer.



Connecting with USB cable





Other features


- Compatible with a wide range of arm circumference (from 17cm to 42cm)
- Compatible with both left and right arm measurements



For right arm

For left arm

Various print-out formats



Commentary sample>
selection mode, the information is printed as follows. troke, and in general do not live as long as people with normal blood pressure. Do not ignore this risk, and monitor your blood pressure regular iv. Regular measurement provides an important means of allerting you to changes in your blood pressure.

Printed measurement results can show not only blood pressure and pulse rate but also knowledge on hypertension.

"name column",
"MAP print", and
"Print QR" are set to
"ON" in the function
selection mode,
the information is

When "Print ID,



Specification

Main unit			
Sales name	OMRON Automatic Blood Pressure Monitor HBP-9030		
General Product Description	Blood pressure monitor for clinical use		
Model	HBP-9030		
Protection Class	Class II		
Operational temperature and humidity	Temperature range : 5 to 40°C Humidity range : 15 to 85%RH (not condensed) Atmospheric pressure : 700 to 1060hPa		
Storage and transportation	Temperature range : -20 to 60°C Humidity range : 10 to 95%RH (not condensed) Atmospheric pressure : 700 to 1060hPa		
Dimensions	Approximately 420 (I) × 460 (w) × 270 (h) mm		
Weight	Approximately 5.5 kg		
Display method	Digital Display		
Accessories included	Power cord, printer paper (1 roll), instruction manual, armrest, AC adapter		

AC adapter

Model	HBP-ACCA-903		
Power	Input Voltage range: AC100 V to 240 V		
	Frequency:	50/60 Hz	

Non-invasive Blood Pressure Measurement Function

Measurement method	Oscillometric method	
NIBP measurement range	Pressure: 0 to 299 mmHg Pulse: 40 to 180 /min	
Accuracy	Pressure: ±3 mmHg Pulse: ±5% of display reading	
Inflation	Automatic inflation	
Measurable arm circumference	17 to 42 cm	

Communication method

USB	USB 2.0	<recommended cable="" usb=""></recommended>
OOD	00D 2.0	
		Comply with USB 2.0, with ferrite core
		at both ends. less than 2.0 m

Separately Sold Items



Armrest Model : HBP-ARM-903

HBP-9030 AC adapter Model: HBP-ACCA-903



Power cord for HBP-9030 AC adapter Model: HBP-CRD-903AP



Model: HBP-9021-CUCO

Cuff cover



Arm Unit with a cuff cover Model: HBP-ARU-903



Printer paper (10 rolls per box) Model: HBP-PAPER-NO16

Distributed by:



Asia Pacific HQ OMRON HEALTHCARE SINGAPORE PTE LTD. 438A Alexandra Road, #05-05/08 Alexandra Technopark, Singapore 119967 www.omronhealthcare-ap.com

Production Facility OMRON DALIAN CO., LTD.

No. 28, Dong Bei Er Street, Economic and Technical Development Zone, Dalian 116600, China

Specifications are subject to change without notice.

• The clinical evaluation of brachial blood pressure measurement was performed in BP mode following ISO 81060-2:2013.

This device is conformed to IEC60601-1-2:2014.